


Digalactosyldiacylglycerol

Chemical Properties

CAS No. :	145033-48-3
Formula:	C49H88O15
Molecular Weight:	917.228
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Digalactosyldiacylglycerol is a natural product for research related to life sciences. The catalog number is TN5814 and the CAS number is 145033-48-3.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0902 mL	5.4512 mL	10.9024 mL
5 mM	0.218 mL	1.0902 mL	2.1805 mL
10 mM	0.109 mL	0.5451 mL	1.0902 mL
50 mM	0.0218 mL	0.109 mL	0.218 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481