


## 1-[2-(2-Pyridyl)ethyl]piperidine

### Chemical Properties

CAS No. : 5452-83-5  
 Formula: C<sub>12</sub>H<sub>18</sub>N<sub>2</sub>  
 Molecular Weight: 190.29  
 Store under nitrogen  
 Storage: Pure form: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



### Biological Description

Description	1-[2-(2-Pyridyl)ethyl]piperidine ,with CAS No. 5452-83-5, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 1-[2-(2-Pyridyl)ethyl]piperidine provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.2551 mL	26.2757 mL	52.5514 mL
5 mM	1.051 mL	5.2551 mL	10.5103 mL
10 mM	0.5255 mL	2.6276 mL	5.2551 mL
50 mM	0.1051 mL	0.5255 mL	1.051 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481