

## 5-[[[(2R,3S,4S,5R,6R)-6-(benzyloxy)-trihydroxyoxan-2-yl]methoxy]-3-hydroxy-methyl-oxopentanoic acid


### Chemical Properties

CAS No. :

Formula: C<sub>19</sub>H<sub>26</sub>O<sub>10</sub>

Molecular Weight: 414.41

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.



### Biological Description

Description	5-[[[(2R,3S,4S,5R,6R)-6-(benzyloxy)-3,4,5-trihydroxyoxan-2-yl]methoxy]-3-hydroxy-3-methyl-5-oxopentanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4131 mL	12.0653 mL	24.1307 mL
5 mM	0.4826 mL	2.4131 mL	4.8261 mL
10 mM	0.2413 mL	1.2065 mL	2.4131 mL
50 mM	0.0483 mL	0.2413 mL	0.4826 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481