


3'-beta-C-Ethynyluridine

Chemical Properties

CAS No. :	180300-49-6
Formula:	C ₁₁ H ₁₂ N ₂ O ₆
Molecular Weight:	268.22
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	3'-beta-C-Ethynyluridine is a Nucleoside Derivative - 3'-Modified nucleoside.
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.7283 mL	18.6414 mL	37.2828 mL
5 mM	0.7457 mL	3.7283 mL	7.4566 mL
10 mM	0.3728 mL	1.8641 mL	3.7283 mL
50 mM	0.0746 mL	0.3728 mL	0.7457 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481