

## Tropanyl trans-cinnamate

## Chemical Properties


CAS No. :

Formula: C<sub>17</sub>H<sub>21</sub>NO<sub>2</sub>

Molecular Weight: 271.36

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Tropanyl trans-cinnamate is a useful organic compound for research related to life sciences and the catalog number is T124966.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6851 mL	18.4257 mL	36.8514 mL
5 mM	0.737 mL	3.6851 mL	7.3703 mL
10 mM	0.3685 mL	1.8426 mL	3.6851 mL
50 mM	0.0737 mL	0.3685 mL	0.737 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481