

29 kDa outer membrane Protein, Acinetobacter baumannii, Recombinant (His & Myc & SUMO)

General Information

Synonyms: 29 kDa outer membrane protein;29 kDa OMP
Protein Construction: 1-24 aa
Species: Acinetobacter baumannii
Expression Host: E. coli
Accession: P83446
Molecular Weight: 18.7 kDa (predicted)
AA Sequence: DEAVVHDSYAFDKNQLIPVGARAE

QC Testing

Biological Activity: Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity: > 85% as determined by SDS-PAGE.
Endotoxin: < 1.0 EU/μg of the protein as determined by the LAL method.
Formulation: Tris-based buffer, 50% glycerol

Preparation and Storage

Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

May be involved in transporting molecules across the outer membrane.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481