

Bovine viral diarrhea virus (strain SD-1) Genome polyprotein (E. coli, His)

General Information

Synonyms:	Genome polyprotein
Protein Construction:	1-168 aa
Species:	BVDV
Expression Host:	E. coli
Accession:	Q01499
Molecular Weight:	25.0 kDa (predicted)
AA Sequence:	MELITNELLYKTYKQKPVGVVEEPVYDQAGNPLFGERGAIHPQSTLKLPHKRGERNVPTSLASLPKRGDCRSGN SKGPVSGIYLKPGPLFYQDYKGPVYHRAPLELFEEGSMCETT KRIGRVTGSDGKLYHIYICIDGCVK SATRSHQ RVLRWVHNRLDCPLWVTSC

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 85% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	Tris-based buffer, 50% glycerol

Preparation and Storage

Reconstitution:
A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Leader cysteine autoprotease that cleaves itself from the nascent polyprotein during translation of the viral mRNA. Once released, plays a role in the inhibition of host innate immune response by interacting with host IRF3 and inducing its proteasomal degradation.; Packages viral RNA to form a viral nucleocapsid and thereby protects viral RNA. Plays also a role in transcription regulation. Protects the incoming virus against IFN-induced effectors.; Initial binding to target cell probably involves interaction of E(rns) with glycosaminoglycans.; E1 and/or E2 are probably

A DRUG SCREENING EXPERT

responsible of cell attachment with CD46 and subsequent fusion after internalization of the virion by endocytosis.; E1 and/or E2 are probably responsible of cell attachment with CD46 and subsequent fusion after internalization of the virion by endocytosis.; Plays an essential role in the virus replication cycle by acting as a viroporin. Forms ion conductive pores, which alters the cell permeability allowing the transport of ions and other small molecules. Forms a leader sequence to properly orient NS2 in the membrane.; Uncleaved NS2-3 is required for production of infectious virus.; Plays a role in the regulation of viral RNA replication.; Multifunctional protein that contains an N-terminal protease and a C-terminal helicase, playing essential roles in viral polyprotein processing and viral genome replication. The chymotrypsin-like serine protease activity utilizes NS4A as an essential cofactor and catalyzes the cleavage of the polyprotein leading to the release of NS4A, NS4B, NS5A, and NS5B. Interacts with NS5B to enhance RNA-dependent RNA polymerase activity.; Acts as a cofactor for the NS3 protease activity.; Induces a specific membrane alteration that serves as a scaffold for the virus replication complex.; Replicates the viral (+) and (-) genome.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481