

AlkB Protein, Caulobacter vibrioides, Recombinant (His & Myc)

General Information

Synonyms:	alkB;DNA oxidative demethylase AlkB;Alpha-ketoglutarate-dependent dioxygenase AlkB homolog
Protein Construction:	1-220 aa
Species:	Caulobacter vibrioides
Expression Host:	E. coli
Accession:	B8GWW6
Molecular Weight:	31.2 kDa (predicted)
AA Sequence:	MAVVRRAVAARGLQMIKPLTVVPGFDVWPGLLDISAQRALVEAVLAGAEQAPFSNYRTAYGKPM SVAMTA LGSLGWTS DARGYRYVDRHPETGRPWPDMPPALLDLWTVLGDPE TTPD SCLVNLYRDGARMGLHQDRDEA DPRFPVLSISLGDTAVFRIGGVNRKDPTRSLRLASGDVCRLLGPARLAFHGVDRILPGSSSLVPGGGGRINLTLRR ARTA

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 85% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	Tris-based buffer, 50% glycerol

Preparation and Storage

Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Dioxygenase that repairs alkylated DNA and RNA containing 3-methylcytosine or 1-methyladenine by oxidative demethylation. Has highest activity towards 3-methylcytosine. Has lower activity towards alkylated DNA containing ethenoadenine, and no detectable activity towards 1-methylguanine or 3-methylthymine. Accepts

A DRUG SCREENING EXPERT

double-stranded and single-stranded substrates. Requires molecular oxygen, alpha-ketoglutarate and iron. Provides extensive resistance to alkylating agents such as MMS and DMS (SN2 agents), but not to MMNG and MNU (SN1 agents).

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481