

MOMP Protein, Chlamydomophila psittaci, Recombinant (His)

General Information

Synonyms: Major outer membrane porin;omp1;MOMP;ompA

Protein Construction: 23-392 aa

Species: Chlamydomophila psittaci

Expression Host: E. coli

Accession: Q00087

Molecular Weight: 45.4 kDa (predicted)

AA Sequence:

```
LPVGNPAEPSLLIDGTMWEGASGDPCDPCATWCDAISIRAGFYGDYVFDRILKVDVNKTISGMAAAPTAASGT
ASNTTVAADRSNFAYGKHLQDAEWCTNAAYLALNIWDRFDVFCFLGASNGYFKASSDAFNLVGLIGLAGTD
FANQRPNVEISQGIVELYTDTAFSWSVGARGALWECGCATLGAEFQYAQSNPKIEMLVNTSSPAQFMIHKPR
GYKGTAAANFPLPVAAGTATATDTKSATVKYHEWQVGLALSRYLNLVPIYGVNWSRATFDADTIRIAQPKLA
SAILNLTWNPTLLGVATLDTSNKYADFMQIVSMQINKMSRKACGIAVGATLIDADKWSITGEARLIDERAA
HINAQFRF
```

QC Testing

Biological Activity: Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.

Purity: > 85% as determined by SDS-PAGE.

Endotoxin: < 1.0 EU/μg of the protein as determined by the LAL method.

Formulation: Tris-based buffer, 50% glycerol

Preparation and Storage

Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months.

Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

In elementary bodies (EBs, the infectious stage, which is able to survive outside the host cell) provides the structural integrity of the outer envelope through disulfide cross-links with the small cysteine-rich protein and the

A DRUG SCREENING EXPERT

large cysteine-rich periplasmic protein. It has been described in publications as the Sarkosyl-insoluble COMC (Chlamydia outer membrane complex), and serves as the functional equivalent of peptidoglycan.; Permits diffusion of specific solutes through the outer membrane.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481