

## SVMP Protein, Crotalus adamanteus, Recombinant (His)

### General Information

Synonyms:	Proteinase II;Adamalysin II;Snake venom metalloproteinase adamalysin-2;SVMP
Protein Construction:	1-203 aa
Species:	Crotalus adamanteus
Expression Host:	E. coli
Accession:	P34179
Molecular Weight:	27.1 kDa (predicted)
AA Sequence:	QQNLPQRYIELVVVADRRVFMKYNSDLNIIIRTRVHEIVNIINGFYRSLNIDVSLVNLEIWSGQDPLTIQSSSSNTL NSEGLWREKVLLNKKKKDNAQLLTAIEFKCETLGKAYLNSMCNPRSSVGIVKDHSPINLLVAVTMAHELGHNL GMEHDGKDCLRGASLCIMRPGLTPGRSYEFSDDSMGYYQKFLNQYKPCILNKP

### QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 90% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	Tris-based buffer, 50% glycerol

### Preparation and Storage

Reconstitution:	A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.
Stability & Storage:	Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots. <small>Actual storage temperature shall be subject to the COA.</small>

### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

### Protein Background

Has no significant hemorrhagic activity, but inactivates serpins by limited proteolysis of their reactive-site loops.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481