

## FtsZ Protein, E. coli O157:H7, Recombinant (His & Myc)

### General Information

Synonyms: ftsZ; Cell division protein FtsZ

Protein Construction: 1-383 aa

Species: E. coli

Expression Host: E. coli

Accession: P0A9A8

Molecular Weight: 47.8 kDa (predicted)

AA Sequence:

```
MFEPMELTNDAVIKVIGVGGGGNAVEHMRERIEGVEFFAVNTDAQALRKTAVGQTIQIGSGITKGLGAGA  
NPEVGRNAADEDRLDALRAALEGADMVFIAAGMGGGTGTGAAPVVAEVAKDLGILTAVVTKPFNFEGKKRM  
AFAEQGITELSKHVDSLITIPNDKLLKVLGRGISLLDAFGAANDVLKGAVQGI AELITRPGLMNVDFADVRTVM  
SEMGYAMMGSVSGEDRAEEAAEMAISPLLEDIDLSGARGVLVNITAGFDLRLDEFETVGN TIRAFASDNA  
TVVIGTSLDPDMNDELRVTVVATGIGMDKRPEITLVTKQVQPVMDRYQQHGMAPLTQE QKPVAKVVND  
NAPQTAKEPDYLDIPAFLRKQAD
```

### QC Testing

Biological Activity: Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.

Purity: > 85% as determined by SDS-PAGE.

Endotoxin: < 1.0 EU/μg of the protein as determined by the LAL method.

Formulation: If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.

### Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in sterile deionized water. The product concentration should not be less than 100 μg/mL. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

*Actual storage temperature shall be subject to the COA.*

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

**Protein Background**

Essential cell division protein that forms a contractile ring structure (Z ring) at the future cell division site. The regulation of the ring assembly controls the timing and the location of cell division. One of the functions of the FtsZ ring is to recruit other cell division proteins to the septum to produce a new cell wall between the dividing cells. Binds GTP and shows GTPase activity.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481