

MscL Protein, E. coli O157:H7, Recombinant (B2M & His)

General Information

Synonyms:	Large-conductance mechanosensitive channel;mscL
Protein Construction:	1-136 aa
Species:	E. coli
Expression Host:	E. coli
Accession:	P0A743
Molecular Weight:	29.0 kDa (predicted)
AA Sequence:	MSIIKEFREFAMRGNVVDLAVGVIIIGAAFGKIVSSLVADIIMPPLGLLIGGIDFKQFAVTLRDAQGDIPAVVMHYGVFIQNVDFDLIVAFMAIKLINKLNRRKKEEPAAAPAPTKEEVLLTEIRDLLKEQNNRS

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 85% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	Tris-based buffer, 50% glycerol

Preparation and Storage

Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Channel that opens in response to stretch forces in the membrane lipid bilayer. May participate in the regulation of osmotic pressure changes within the cell.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481