

YcbX Protein, E. coli, Recombinant (GST)

General Information

Synonyms: ycbX;Uncharacterized protein YcbX

Protein Construction: 1-369 aa

Species: E. coli

Expression Host: E. coli

Accession: P75863

Molecular Weight: 67.6 kDa (predicted)

AA Sequence:

MATLIRLFHPVKSMRGI GLTHALADV SGLAFDRIFMITEPDGTFITARQFPQMVRFTSPSPVHDGLHLTAPDGSS
AYVRFADFATQDAPTEVWGTHFTARIAPDAINKWLSGFFSREVQLRWVGPQMTRRVK RHNTVPLSFADGYP
YLLANEASLRDLQQRCPASVKMEQFRPNLVVSGASAW EEDRWK VIRIGDVVFDVVKPCSRCIFTTVSPEKGQ
KHPAGEPLKTLQSFRTAQDNGDVDFGQNL IARNSGVIRVGDVEILATAPAKIYGAAAADDTANITQQPDAN
VDIDWQGQAFRGN NQQVLLEQLENQ GIRIPYSCRAGICGSCR VQ LLEGEVTPLKKSAMGDDG TILCCSCVPKT
ALKLAR

QC Testing

Biological Activity: Measured by its binding ability in a functional ELISA. Immobilized IdC at 2 µg/mL can bind human ycbX, the EC50 of human ycbX protein is 40.54-47.97 µg/mL.

Purity: > 85% as determined by SDS-PAGE.

Endotoxin: < 1.0 EU/µg of the protein as determined by the LAL method.

Formulation: Tris-based buffer, 50% glycerol

Preparation and Storage

Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months.

Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

N/A

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481