

ABAT Protein, Human, Recombinant (His & Myc & SUMO)

General Information

Synonyms:	GABAT;4-aminobutyrate aminotransferase, mitochondrial;GABA aminotransferase (GABA-AT);(S)-3-amino-2-methylpropionate transaminase;ABAT;L-AIBAT;Gamma-amino-N-butyrate transaminase (GABA transaminase;GABA-T);mitochondrial
Protein Construction:	29-500 aa
Species:	Human
Expression Host:	E. coli
Accession:	P80404
Molecular Weight:	73.3 kDa (predicted)
AA Sequence:	SQAAAKVDVEFDYDGPLMKTEVPGPRSQELMKQLNIIQNAEAVHFFCNYYESRGNLVDVDGNRMLDLYSQ ISSVPIGYSHPLLKLIQPPQNASMFVNRPALGILPPENFVEKLRQSLLSVAPKGMSQLITMACGSCSNENALK TIFMWYRSKERGQRGFSQEELETMINQAPGCPDYSILSFMGAFHGRTMGCLATTHSKAIHKIDIPFDWPIAP FPRLKYPLEEFVKENQQEEARCLEEVEDLIVKYRKKKKTVAGIIVEPIQSEGDNHASDDFFRKLRLDIARKHGCA FLVDEVQTGGGCTGKFWAHEHWGLDDPADVMTFSKKMMTGFFHKEEFRPNAPYRIFNTWLGDPKSNLLL AEVINIIRKEDLLNNAAHAGKALLTGLLDLQARYPQFISVRVGRGTFCSDTPDDSI RNKLIL IARNKGVVLGGC GDKSIRFRPTLVFRDHHHLFLNIFSDILADFK

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 90% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	Tris-based buffer, 50% glycerol

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 μg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Catalyzes the conversion of gamma-aminobutyrate and L-beta-aminoisobutyrate to succinate semialdehyde and methylmalonate semialdehyde, respectively. Can also convert delta-aminovalerate and beta-alanine.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481