

## CIRBP Protein, Human, Recombinant (GST)

### General Information

Synonyms:	A18HNRNP;Glycine-rich RNA-binding protein CIRP;A18 hnRNP;CIRP;CIRBP;Cold-inducible RNA-binding protein
Protein Construction:	1-172 aa
Species:	Human
Expression Host:	E. coli
Accession:	Q14011
Molecular Weight:	45.6 kDa (predicted)
AA Sequence:	MASDEGKLFVGGLSFDTNEQSLEQVFSKYGQISEVVVVVKDRETQRSRGGFVTFENIDDAKDAMMAMNGKS VDGRQIRVDQAGKSSDNRSRGYRGGGAGGRGFRGGRGRGFRGGGDRGYGGNRFESRGGYGGSRDY YSSRSQSGGYSRSGGYSYRDSYDSYATHNE

### QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 90% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	Tris-based buffer, 50% glycerol

### Preparation and Storage

#### Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

#### Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

### Protein Background

Cold-inducible mRNA binding protein that plays a protective role in the genotoxic stress response by stabilizing transcripts of genes involved in cell survival. Acts as a translational activator. Seems to play an essential role in cold-induced suppression of cell proliferation. Binds specifically to the 3'-untranslated regions (3'-UTRs) of stress-responsive transcripts RPA2 and TXN. Acts as a translational repressor. Promotes assembly of stress granules

(SGs), when overexpressed.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286    E\_mail: info@targetmol.com    Address: 34 Washington Street, Wellesley Hills, MA 02481