

DUSP26 Protein, Human, Recombinant (His & SUMO)

General Information

Synonyms:	Low-molecular-mass dual-specificity phosphatase 4 (DSP-4;LDP-4);MKP8;Novel amplified gene in thyroid anaplastic cancer;LDP4;DUSP24;Dual specificity protein phosphatase 26; NATA1;DUSP26;Dual specificity phosphatase SKRP3;SKRP3;Mitogen-activated protein kinase phosphatase 8 (MAP kinase phosphatase 8;MKP-8)
Protein Construction:	1-211 aa
Species:	Human
Expression Host:	E. coli
Accession:	Q9BV47
Molecular Weight:	39.9 kDa (predicted)
AA Sequence:	MCPGNLWASMTFMARFSRSSRSPVTRTRGTLEEMPTVQHPFLNVFELERLLYTGKTACNHADEVWPGLYL GDQDMANNRRELRRLLGITHVLNASHSRWRGTPEAYEGLGIRYLGVEAHDSPAFDMSIHFTAADFIHRALSQ PGGKILVHCAVGVSRSATLVLAYLMLYHHLTLVEAIKKVKDHRGIIPNRGFLRQLLALDRRLRQGLEA

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 90% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	Tris-based buffer, 50% glycerol

Preparation and Storage

Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Inactivates MAPK1 and MAPK3 which leads to dephosphorylation of heat shock factor protein 4 and a reduction in its DNA-binding activity. Inhibits MAP kinase p38 by dephosphorylating it and inhibits p38-mediated apoptosis in

anaplastic thyroid cancer cells. Can also induce activation of MAP kinase p38 and c-Jun N-terminal kinase (JNK).

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481