

## EIF4A2 Protein, Human, Recombinant (GST)

### General Information

Synonyms:	eIF4A-II; eIF-4A-II; DDX2B; EIF4F; EIF4A2; ATP-dependent RNA helicase eIF4A-2; Eukaryotic initiation factor 4A-II
Protein Construction:	1-407 aa
Species:	Human
Expression Host:	E. coli
Accession:	Q14240
Molecular Weight:	73.4 kDa (predicted)
AA Sequence:	<p>MSGGSADYNREHGGPEGMDPDGVIESNWNEIVDNFDDMNLKESLLRGIYAYGF EKPSAIQQR AII PCIKGYD VIAQAQSGTGKTATFAISILQQLIEFKETQALVLA PTRELAQQIQVILALGDYMGATCHACIGGTNVRNEMQ KLQAEAPHIVGTPGRVFDMLNRRYLSPKWIKMFVLDEADEMLSRGFKDQIYEIFQKLNTSIQVVLLSATMPT DVLEVTKKFM RDP IRLVKKEELTLEGIKQFYINVEREEWKLDTLCDLYETLTITQAVIFLNTRRKVDWLTEKMHA RDFTVSALHGDM DQKERDVIMREFRSGSSRVLITDLLARGIDVQVSLVINYLPTNRENYIHRIGRGGFRGFR KGVAINFVTEEDKRILRDIETFYNTTVEEMPMNVADLI</p>

### QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 90% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	Lyophilized from PBS, 6% Trehalose, pH 7.4

### Preparation and Storage

#### Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

#### Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481