

FGFR1 Protein, Human, Recombinant (E. coli, His)

General Information

Synonyms:	Proto-oncogene c-Fgr;Fms-like tyrosine kinase 2 (FLT-2);FLG;BFGFR;Basic fibroblast growth factor receptor 1 (BFGFR;bFGF-R-1);CD331;N-sam;FGFBR;FLT2;FGFR-1;HBGFR;CEK;FGFR1; Fibroblast growth factor receptor 1
Protein Construction:	22-376 aa
Species:	Human
Expression Host:	E. coli
Accession:	P11362
Molecular Weight:	43.4 kDa (predicted)
AA Sequence:	RPSPTLPEQAQPWGAPVEVESFLVHPGDLLQLRCRLRDDVQINWLRDGVQLAESNRTRITGEEVEVQDSVP ADSGLYACVTSSPSGSDTTYFSVNVSDALPSEDDDDDDSSSEEKETDNTKPNRMPVAPYWTSPEKMEKK LHAVPAAKTVKFKCPSSGTPNPTLRWLKNGKEFKPDHRIGGYKVRYATWSIIMDSVVPSPDKGNVTCIVENEYG SINHLYQLDVVERSPHRPILQAGLPANKTVALGSNVEFMCKVYSDPQPHIQWLKHIEVNGSKIGPDNLPYVQI LKTAGVNTTDKEMEVLHLRNVSFEDAGEYTCLAGNSIGLSHSAWLTVLEALEERPAVMTSPLYLE

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 85% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	Tris-based buffer, 50% glycerol

Preparation and Storage

Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Tyrosine-protein kinase that acts as cell-surface receptor for fibroblast growth factors and plays an essential role

in the regulation of embryonic development, cell proliferation, differentiation and migration. Required for normal mesoderm patterning and correct axial organization during embryonic development, normal skeletogenesis and normal development of the gonadotropin-releasing hormone (GnRH) neuronal system. Phosphorylates PLCG1, FRS2, GAB1 and SHB. Ligand binding leads to the activation of several signaling cascades. Activation of PLCG1 leads to the production of the cellular signaling molecules diacylglycerol and inositol 1,4,5-trisphosphate. Phosphorylation of FRS2 triggers recruitment of GRB2, GAB1, PIK3R1 and SOS1, and mediates activation of RAS, MAPK1/ERK2, MAPK3/ERK1 and the MAP kinase signaling pathway, as well as of the AKT1 signaling pathway. Promotes phosphorylation of SHC1, STAT1 and PTPN11/SHP2. In the nucleus, enhances RPS6KA1 and CREB1 activity and contributes to the regulation of transcription. FGFR1 signaling is down-regulated by IL17RD/SEF, and by FGFR1 ubiquitination, internalization and degradation.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481