

GAL3ST3 Protein, Human, Recombinant (His & Myc)

General Information

Synonyms:	Gal3ST-3;Galactose-3-O-sulfotransferase 3;Beta-galactose-3-O-sulfotransferase 3 (Gal3ST3);Gal-beta-1, 3-GalNAc 3'-sulfotransferase 3;GAL3ST3
Protein Construction:	41-431 aa
Species:	Human
Expression Host:	E. coli
Accession:	Q96A11
Molecular Weight:	51.9 kDa (predicted)
AA Sequence:	SWYPKLFPLSCPPLRNSPPRPKHMTVAFLKTHKTAGTTVQNILFRFAERHNLTVALPHPSCEHQFCYPRNFSAHFVHPATRPPHVLASHLRFDRALERLMPPSTVYVITLREPAAMFESLFSYYNQYCPAFRRVNPASLEAFLRAP EAYYRAGEHFAMFAHNTLAYDLGGDNERSPRDDAAYLAGLIRQVEEVFSLVMIAEYFDESLVLLRLLAWDL DDVLYAKLNARAASSRLAAIPAALARAARTWNALDAGLYDHFNATFWRHVARAGRACVEREARELREARQ RLLRRCFGDEPLLRPAAQIRTKQLQPWQPSRKVDIMGYDLPGGGAGPATEACLKLMPEVQYSNYLLRKQKR RGGARARPEPVLNPPPRPIRVLPRGPQGP

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 85% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in sterile deionized water. The product concentration should not be less than 100 μg/mL. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Transfers a sulfate to position 3 of non-reducing beta-galactosyl residues in N-glycans and core2-branched O-glycans. Has high activity towards Gal-beta-1,4-GlcNAc, Gal-beta-1,4(Fuc-alpha-1,3)GlcNAc and lower activity towards Gal-beta-1,3(Fuc-alpha-1,4)GlcNAc.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481