

Human herpesvirus 6A (HHV-6 variant A) (strain Uganda-1102) DNA polymerase processivity

General Information

Synonyms:	Polymerase accessory protein (PAP);Phosphoprotein P41 (PP41);U27;DNA polymerase processivity factor;EPLF1
Protein Construction:	1-393 aa
Species:	HHV-6A
Expression Host:	P. pastoris (Yeast)
Accession:	P52439
Molecular Weight:	46.8 kDa (predicted)
AA Sequence:	MCWSFHLLFFKAHKARVGGARTSFLTEMERGSRDHHRDHRDHREHRETREPPTLAFHMKSWKTINKSLKAFK LLKENTTVTFTPQPSIIISAKNHLVQKLTIQAECLFLSDTDRFLTKTINNHIFESFMNIISNPEVTKMYIQHDS DLYTRVLVTASDTCTQASVPCVHGQEVVRDTGRSPLRIDLDHSTVSDVLKWLSPVTKTKRSGKSDALMAHIIV QVNPPTIKFVTEMNELEFSNSNKVIFYDVKNMRFNLSAKNLQQALSMCAVIKTSCSLRTVAAKDCKLILTSKST LLTVEAFLTQEQLKEESRFERMGKQDDGKGDGRSHKNDGGSALASKQEMQYKITNYMVPKNGTAGSSLFNE KEDSESDDSMHFDYSSNPNPKRQRCVV

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 90% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	Tris-based buffer, 50% glycerol

Preparation and Storage

Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Accessory subunit of the DNA polymerase that acts to increase the processivity of polymerization.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481