

## Human herpesvirus 6A (HHV-6 variant A) (strain Uganda-1102) Envelope glycoprotein B (His &

### General Information

Synonyms:	gB;Envelope glycoprotein B
Protein Construction:	23-188 aa
Species:	HHV-6A
Expression Host:	E. coli
Accession:	P28864
Molecular Weight:	26.8 kDa (predicted)
AA Sequence:	DPDHYIRAGYNHKYPFRICSIAKGTDLMRFRDRDISCSPYKSNAKMSEGFFIYKTNIEYTFPVRTYKKELTFQSS YRDVGVVYFLDRTVMGLAMPVYEANLVNSHAQCYSAVAMKRPDGTVFSAFHEDNNKNNTLNLPLNFKSIT NKRFITTKEPYFARGPLWL

### QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 85% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.

### Preparation and Storage

#### Reconstitution:

Reconstitute the lyophilized protein in sterile deionized water. The product concentration should not be less than 100 μg/mL. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

#### Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

### Protein Background

Envelope glycoprotein that forms spikes at the surface of virion envelope. Essential for the initial attachment to

heparan sulfate moieties of the host cell surface proteoglycans. Involved in fusion of viral and cellular membranes leading to virus entry into the host cell. Following initial binding to its host receptors, membrane fusion is mediated by the fusion machinery composed at least of gB and the heterodimer gH/gL. May be involved in the fusion between the virion envelope and the outer nuclear membrane during virion egress.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481