

HABP2 Protein, Human, Recombinant (His & SUMO)

General Information

Synonyms:	Hyaluronan-binding protein 2;PHBP;Factor VII-activating protease;Plasma hyaluronan-binding protein;Hepatocyte growth factor activator-like protein;Factor seven-activating protease (FSAP);HABP2;HGFAL
Protein Construction:	314-560 aa
Species:	Human
Expression Host:	E. coli
Accession:	Q14520
Molecular Weight:	43.3 kDa (predicted)
AA Sequence:	IYGGFKSTAGKHPWQASLQSSLPLTISMPQGHFCGGALIHPCWVLTAAHCTDIKTRHLKVVLDQDLKKEEFH EQSFRVEKIFKYSHYNERDEIPHNDIALLLKLPVDGHCALSKYVKTVCLPDGSPSGSECHISGWGVTTETGKG SRQLLDKVKLIANTLCSNRQLYDHMIDDSMICAGNLQKPGQDTCQGDSGGPLTCEKDGTYVYVYGVISWGLE CGKRPGVYTQVTKFLNWIKATIKSESGF

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 90% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	Tris-based buffer, 50% glycerol

Preparation and Storage

Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Cleaves the alpha-chain at multiple sites and the beta-chain between 'Lys-53' and 'Lys-54' but not the gamma-chain of fibrinogen and therefore does not initiate the formation of the fibrin clot and does not cause the

A DRUG SCREENING EXPERT

fibrinolysis directly. It does not cleave (activate) prothrombin and plasminogen but converts the inactive single chain urinary plasminogen activator (pro-urokinase) to the active two chain form. Activates coagulation factor VII. May function as a tumor suppressor negatively regulating cell proliferation and cell migration.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481