

## KRT4 Protein, Human, Recombinant (His & SUMO)

### General Information

Synonyms:	Cytokeratin-4 (CK-4);CYK4;KRT4;Keratin, type II cytoskeletal 4;Type-II keratin Kb4;type II cytoskeletal 4;Keratin-4 (K4)
Protein Construction:	1-534 aa
Species:	Human
Expression Host:	E. coli
Accession:	P19013
Molecular Weight:	73.3 kDa (predicted)
AA Sequence:	MIARQQCVRGGPRGFSCGSAIVGGGKRGAFSSVSMMSGGAGRCSSGGFGSRSLYNLRGNKSISMSVAGSRQGACFGGAGGFGTGGFGAGGFGAGFGTGGFGGGFGGSGFSGKGGPGFPVCPAGGIQEVNTINQSLTPLHVEIDPEIQKVRTEEREQIKLLNNKFASFIDKVQFLEQQNKVLETKWNLLQQQTTTTSSKNLEPLFETYLSVLRKQLDTLGN DKGRLQSELKTMQDSVEDFKTKYEEEEINKRTAAENDFVVLKKDVAAYLNKVELEAKVDSLNDINFLKVLVD AELSQMQTHVSDTSVVLSDMNNRNLDLDSIIAEVRAQYEEIAQRSKAEAEALYQTKVQQQLQISVDQHGDNLK NTKSEIAELNRMQRLRAEIEINIKKQCQTLQVSVADAEQRGENALKDAHSKRVELEAALQAKEELARMLREY QELMSVKLALDIEIATYRKLLEGEYRMSGECQSAVSISVSVSGSTSTGGISGGLGSGSGFGLSSGFGSGSGSGFG FGGSVSGSSSSKIISTTTLNKR

### QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 90% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	Tris-based buffer, 50% glycerol

### Preparation and Storage

#### Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

#### Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

*Actual storage temperature shall be subject to the COA.*

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

N/A

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481