

CT83 Protein-VLP, Human, Recombinant (His)

General Information

Synonyms:	Cancer/testis antigen 83;CT83;KK-LC-1;CXorf61;Kita-kyushu lung cancer antigen 1;KKLC1
Protein Construction:	1-113 aa
Species:	Human
Expression Host:	HEK293 Cells
Accession:	Q5H943
Molecular Weight:	14.0 kDa (predicted)
AA Sequence:	MNFYLLASSILCALIVFWKYRRFQRNTGEMSSNSTALALVRPSSSGLINSNTDNNLAVYDLSDILNNFPHSI ARQKRILVNLSMVENKLVELEHTLLSKGFRGASPHRKST

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	PBS, pH 7.4, 6% trehalose

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in sterile deionized water. The product concentration should not be less than 100 μg/mL. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

N/A

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481