

PMEL Protein, Human, Recombinant (His & SUMO)

General Information

Synonyms:	Silver locus protein homolog;Premelanosome protein;ME20-M (ME20M);Melanoma-associated ME20 antigen;D12S53E;Melanocyte protein PMEL;Melanocyte protein Pmel 17;PMEL;P1;PMEL17;Melanocytes lineage-specific antigen GP100;SILV;P100
Protein Construction:	25-467 aa
Species:	Human
Expression Host:	E. coli
Accession:	P40967
Molecular Weight:	63.0 kDa (predicted)
AA Sequence:	KVPRNQDWLGVSRQLRTKAWNRQLYPEWTEAQLDCWRGGQVSLKVSNDGPTLIGANASFSIALNFPGSQ KVLDPDGQVIWVNNTIINGSQVWGGQPVPYQETDDACIFPDGGPCPSGSWSQKRSFVYVWKTWGQYWQVL GGPVSGLSIGTGRAMLGHTMEVTYHRRGSRYSVPLAHSSSAFTITDQVPFVSVSQRLRALDGGNKHFLRN QPLTFALQLHDPSGYLAEADLSYTWDFGDSSGTLISRALVVTHTYLEPGPVTAQVVLQAAIPLTSCGSSPVPPT TDGHRPTAEAPNTTAGQVPTTEVVGTTGQAPTAEPSGTTSVQVPTTEVISTAPVQMPTAESTGMTPEKVPVS EVMGTTLAEMSTPEATGMTPAEVSIVVLSGTTAAQVTTTEWVETTARELPIPEPEGPDASSIMSTESITGSLGPL LDGTATLRLV

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 90% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	Tris-based buffer, 50% glycerol

Preparation and Storage

Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Plays a central role in the biogenesis of melanosomes. Involved in the maturation of melanosomes from stage I to II. The transition from stage I melanosomes to stage II melanosomes involves an elongation of the vesicle, and the appearance within of distinct fibrillar structures. Release of the soluble form, ME20-S, could protect tumor cells from antibody mediated immunity.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481