

## MYBPC2 Protein, Human, Recombinant (His)

## General Information

Synonyms:	fast-type;C-protein, skeletal muscle fast isoform;MYBPC2;MYBPCF;Myosin-binding protein C, fast-type;Fast MyBP-C
Protein Construction:	739-1141 aa
Species:	Human
Expression Host:	E. coli
Accession:	Q14324
Molecular Weight:	49.5 kDa (predicted)
AA Sequence:	EPLHLIVEDVTDTTTLKWRPPNRIGAGGIDGYLVEYCLEGSEEWVPANTEPVERCGFTVKNLPTGARILFRVVG VNIAGRSEPATLAQPVTIREIAEPPKIRLPRHLRQTYIRKVGQNLNLVVPFQKPRPQVVWTKGGAPLDTSRVH VRTSDFDTVFFVRQAARSDSGEYELSVQIENMKDTATIRIRVVEKAGPPINVMVKEVWGTNALVEWQAPKDD GNSEIMGYFVQKADKKTMEWFNVYERNRHTSCTVSDLIVGNEYFRVYTENICGLSDSPGVSKNTARILKTGIT FKPFYKEHDFRMAPKFLTPLIDRVVAGYSAALNCAVRGHPKPKVWVMKNKMEIREDPKFLITNYQGVLL NIRRSPFDAGTYTCRAVNELGEALAECKLEVRVPQ

## QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 90% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	Tris-based buffer, 50% glycerol

## Preparation and Storage

## Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

## Stability &amp; Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

## Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

## Protein Background

Thick filament-associated protein located in the crossbridge region of vertebrate striated muscle a bands. In vitro

## A DRUG SCREENING EXPERT

---

it binds MHC, F-actin and native thin filaments, and modifies the activity of actin-activated myosin ATPase. It may modulate muscle contraction or may play a more structural role.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481