

## CYPOR Protein, Human, Recombinant (GST)

### General Information

Synonyms: NADPH--cytochrome P450 reductase;CYPOR;CPR;P450R;POR

Protein Construction: 2-671 aa

Species: Human

Expression Host: E. coli

Accession: P16435

Molecular Weight: 102.9 kDa (predicted)

AA Sequence:

```
INMGDSHVDTSSVSEAVAEVSLFSMTDMILFSLIVGLLTYWFLFRKKKEEVPEFTKIQTLTSSVRESSFVEKMK
KTGRNIIVFYGSQTGTAEFANRLSKDAHRYGMRGMSADPEEYDLADLSSLPEIDNALVVFVCMATYGEDPT
DNAQDFYDWLQETDVDLSGVKFAVFGNGKTYEHFNAMGKYVDRLEQLGAQRIFELGLGDDDDGNLEEDFI
TWREQFWLAVCEHFGVEATGEESSIRQYELVVHTDIDAAKVYMGEMGRKLSYENQKPPFDAKNPFLAAVTTN
RKLNQGTERHLMHLELDISDSKIRYESGDHVAVYPANDSALVNQLGKILGADLDVVMNLNDEESNKKHPF
PCPTSYRTALTYLDITNPPRTNVLYELAQYASEPSEQELLRKMASSGEGKELYLSWVVEARRHILAILQDCPS
LRPPIDHLCCELLPRLQARYYSIASSSKVHPNSVHICAVVVEYETKAGRINKGVATNWLRAKEPVGENGGRALVP
MFVRKSFRLPFKATTPVIMVGPVGTGVPFIGFIQERAWLRQQGKEVGETLLYYGCRRSDEDYLYRELAQFH
RDGALTQLNVAFSREQSHKVYVQHLLKQDREHLWKLIEGGAHIYVCGDARNMARDVQNTFYDIVAELGAME
HAQAVDYIKKLMTK
```

### QC Testing

Biological Activity: Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.

Purity: > 90% as determined by SDS-PAGE.

Endotoxin: < 1.0 EU/μg of the protein as determined by the LAL method.

Formulation: Tris-based buffer, 50% glycerol

### Preparation and Storage

Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

**Protein Background**

This enzyme is required for electron transfer from NADP to cytochrome P450 in microsomes. It can also provide electron transfer to heme oxygenase and cytochrome B5.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481