

## Human parainfluenza 4b virus (HPIV-4b) Nucleoprotein/NP Protein (His)

## General Information

Synonyms: NP;Nucleoprotein;Nucleocapsid protein (NP;Protein N)

Protein Construction: 1-551 aa

Species: HPIV-4b

Expression Host: E. coli

Accession: P17241

Molecular Weight: 68.5 kDa (predicted)

AA Sequence: MSSVLAAYEQFLQTTEDRSFGDQQFVQSDTLKAEIPVFLNTNDPQQRFTLMNFCLRLAVSSSAKSAIKQGAL  
LSLLSLQATSMQNHLMIAARAPDAALRIIEVDAIDPQDYTLTINPRSGWDDIKIRAYRALSRLPISLADRTV  
SRDAEHAVCDDMDTYLNRIFSVLIQIWIMVCKCMTAYDQPTGSEERRLAKYKQQGRMLEKYQLQTDARKIIQL  
VIRESMIRQFLVQEMLTADKVGAYTNRYYAMVGDIKYIANVGMSAFFLTKFGLGNRWKPLALAAFGELV  
KLKSFMSLYRRLGDRSRYLALLEPELMEFAPANYPYLLFSYAMGVGSVQDPLIRNYQFGRNFLNTSYFYGVE  
TAMKHQGTVDPKLALLEGITDEDVDIMQSVKHSKAGDDISQAGAFMTMSLSRSFAFINNNTSQDFSGARL  
SNYEQGWSGTNQDETRDVPESTMHRLQNIPTSDHNEMLMPELESDSNPFNRPRFTVRAPLIPEISHQNP  
TTRMNRNINTRDNTRADHQDTNEDRGSNVPDDILGDLN

## QC Testing

Biological Activity: Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.

Purity: > 85% as determined by SDS-PAGE.

Endotoxin: < 1.0 EU/μg of the protein as determined by the LAL method.

Formulation: If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.

## Preparation and Storage

## Reconstitution:

Reconstitute the lyophilized protein in sterile deionized water. The product concentration should not be less than 100 μg/mL. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

## Stability &amp; Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

*Actual storage temperature shall be subject to the COA.*

## Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

**Protein Background**

Encapsidates the genome, protecting it from nucleases. The nucleocapsid (NC) has a helical structure. The encapsidated genomic RNA is termed the NC and serves as template for transcription and replication. During replication, encapsidation by N is coupled to RNA synthesis and all replicative products are resistant to nucleases.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481