

## PARP2 Protein, Human, Recombinant (GST & His)

### General Information

Synonyms:	PARP-2;DNA ADP-ribosyltransferase PARP2;ADPRT2;NAD(+) ADP-ribosyltransferase 2 (ADPRT-2);ADP-ribosyltransferase diphtheria toxin-like 2 (ARTD2);Poly[ADP-ribose] synthase 2 (pADPRT-2);Poly [ADP-ribose] polymerase 2;Protein poly-ADP-ribosyltransferase PARP2; hPARP-2;ADPRTL2;PARP2
Protein Construction:	1-583 aa
Species:	Human
Expression Host:	Baculovirus Insect Cells
Accession:	Q9UGN5
Molecular Weight:	92.8 kDa (predicted)
AA Sequence:	MAARRRRSTGGGRARALNESKRVNNGNTAPEDSSPAKKTRRCQRQESKKMPVAGGKANKDRTEDKQDGM PGRSWASKRVSESVKALLLKGPVDPECTAKVGVKAHVYCEGNDVYDVMLNQTNLQFNNNKYLIQLLEDD AQRNFSVWMRWGRVGMKGQHSVLVACSGNLNKAKEIFQKKFLDKTKNNWEDREKFEKVPKGYDMLQMDY ATNTQDEEETKKEESLKSPLKPESQLDLRVQELIKLICNVQAMEEMMEMKYNTKKAPLGKLTVAQIKAGYQS LKKIEDCIRAGQHGRALMEACNEFYTRIPHDFGLRTPPLIRTQKELSEKIQLLEALGDIEIAIKLVKTELQSPEHPL DQHYRNLHICALRPLDHESYEFKVISQYLQSTHAPTHSDYTMLLDLFEVEKDGEKEAFREDLHNRMLLWHGS RMSNWVGLSHGLRIAPPEAPITGYMFGKGIYFADMSSKSANYCFASRLKNTGLLLLSEVALGQCNELLEANP KAEGLLQGHSTKGLGKMAPSSAHFVTLNGSTVPLGPASDTGILNPDGYTLNNEYIVYNPNQVRMRYLLKV QFNFLQLW

### QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 85% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	Lyophilized from 10 mM Tris-HCL, 1 mM EDTA, 6% Trehalose, pH 8.0

### Preparation and Storage

#### Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

#### Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months.

Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481