

Cereblon/CRBN Protein, Human, Recombinant (Cell-Free, His)

General Information

Synonyms: CRBN;Protein cereblon

Protein Construction: 1-442 aa

Species: Human

Expression Host: E. coli

Accession: Q96SW2

Molecular Weight: 51.2 kDa (predicted)

MAGEGDQQDAAHNMGHNLPLPAESEEEDEMEVEDQDSKEAKKPNIIINFDTSLPTSHTYLGADMEEFHGRT
LHDDDCSQVIPVLPQVMMILIPGQTLPLQLFHPQEVSMVRNLIQKDRTFVLAYSINVQEREAQFGTTAEIYAY
REEQDFGIEIVKKAIGRQRFKVLRLTQSDGIQQAKVQILPECVLPSTMSAVQLESLNKCQIFPSKPVSREDQC
AA Sequence: SYKWWQKYQKRKFHCANLTSWPRWLYSLYDAETLMDRIKKQLREWDENLKDDSLPSNPIDFSYRVAACLPI
DDVLRIQLLKIGSAIQRLRCELDIMNKCTSLCKQCQETEITTKNEIFSLSLCGPMAAYVNPVPHGYVHETLVYKA
CNLNLIGRPSTEHSWFPGYAWTVAQCKICASHIGWKFTATKKDMSPQKFWGLTRSALLPTIPDTEDEISPDKVI
LCL

QC Testing

Biological Activity: Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.

Purity: > 90% as determined by SDS-PAGE.

Endotoxin: < 1.0 EU/μg of the protein as determined by the LAL method.

Formulation: Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Preparation and Storage

Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

A DRUG SCREENING EXPERT

Substrate recognition component of a DCX (DDB1-CUL4-X-box) E3 protein ligase complex that mediates the ubiquitination and subsequent proteasomal degradation of target proteins, such as MEIS2 (Probable). Normal degradation of key regulatory proteins is required for normal limb outgrowth and expression of the fibroblast growth factor FGF8. Maintains presynaptic glutamate release and consequently cognitive functions, such as memory and learning, by negatively regulating large-conductance calcium-activated potassium (BK) channels in excitatory neurons. Likely to function by regulating the assembly and neuronal surface expression of BK channels via its interaction with KCNT1. May also be involved in regulating anxiety-like behaviors via a BK channel-independent mechanism.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481