

RNF207 Protein, Human, Recombinant (His & Myc)

General Information

Synonyms: RNF207;RING finger protein 207

Protein Construction: 1-634 aa

Species: Human

Expression Host: E. coli

Accession: Q6ZRF8

Molecular Weight: 78.3 kDa (predicted)

AA Sequence: MSGAIFGPLEGPSSLDAPSIHPLVCPLCHVQYERPCLLDCFHDFCAGCLRGRATDGR LTCPLCQHQT V LKGPS
GLPPVDRLLQFLVDSSGDGVEAVRCANCDLECSEQDVETTYFCNTCGQPLCARCRDETHRARMFARH DIVAL
GQRSRDVPQKCTLHAEPYLLFSTDKKLLLCIRCFRDMQKESRAHCVDLESAYVQG CERLEQAVLAVKALQTAT
REAIALLQAMVEEV R HSAAEEEDAIHALFGSMQDR LAERKALLLQAVQSQYEEKDKAFKEQLSHLATLLPTLQ
VHLVICSSFLSLANKAEFLDLGYELMERLQGIVTRPHHLRPIQSSKIASDHRAEFARCLEPLLLL G PRRVAAAAS
GANTLAGGLGPKALTGPHCPSPVKGMSGSPVQKPTLHR S I STKVLLAEGENTPFAEHCRHYEDSYRHLQAEM
QSLKDQVQELHRDLTKHHS LIKAEIMGDV L H KSLQLDVQIASEHASLEGMRVVFQEIWEEAYQRVANEQEIYE
AQLHDLLQLRQENAYLTTITKQITPYVRSIAKVKERLEPRFQAPVDEQSESLQNT H D DSRNNAASARNNPGSV
PEKREKTSEPKGNSWAPNGLSEEP L LKNMDHHR SKQKNGGDVPTWREHPT

QC Testing

Biological Activity: Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.

Purity: > 85% as determined by SDS-PAGE.

Endotoxin: < 1.0 EU/μg of the protein as determined by the LAL method.

Formulation: If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in sterile deionized water. The product concentration should not be less than 100 μg/mL. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months.

Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Plays a role in cardiac repolarization possibly by stabilizing membrane expression of the potassium channel KCNH2/HERG, or by assisting its synthesis, folding or export from the endoplasmic reticulum, in a heat shock protein-dependent manner.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481