

SLFN12 Protein, Human, Recombinant (His)

General Information

Synonyms: Schlafen family member 12; Ribonuclease SLFN12; SLFN12

Protein Construction: 1-578 aa

Species: Human

Expression Host: E. coli

Accession: Q8IYM2

Molecular Weight: 71.0 kDa (predicted)

AA Sequence: MNISVDLETNYAELVLDVGRVTLGENSRKKMKDCKLRKKQNESVSRAMCALLNSGGGVKAEIENEDYSYTK
DGIGLDLENSFSNILLFVPEYLDPMQNGNYFLIFVKSWSLNTSGLRITTLSSNLYKRDITSAKVMNATAALEFLK
DMKKTRGRLYLRPELLAKRPCVDIQEENNMKALAGVFFDRTELDRKEKLTFTTESTHVEIKNFSTEKLLQRIKEILP
QYVSFAFANTDGGYLFIGNEDKEIIGFKAEMSDLDDLREIEKSIRKMPVHHFCMEKKKINYSCKFLGVYDKGS
LCGYVCALRVERFCCAFAKEPDSWHVKDNRVMQLTRKEWIQFMVEAEPKFSSSYEEVISQINTSLPAPHSW
PLLEWQRQRHHCPGLSGRITYTPENLCRKLFLQHEGLKQLICEEMDSVRKGSIFSRWSVDLGLQENHKVLC
DALLISQDSPPLYTFHMQDEEFKGYSTQATLTKQKLAIGGYTKKVCVMTKIFYLSPEGMTSCQYDLRSQV
IYPESYFTRRKYLLKALFKALKRLKSLRDQFSFAENLYQIIGIDCFQKNDKMKFKSCRRLT

QC Testing

Biological Activity: Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.

Purity: > 90% as determined by SDS-PAGE.

Endotoxin: < 1.0 EU/μg of the protein as determined by the LAL method.

Formulation: Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Preparation and Storage

Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

N/A

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481