

SLC23A2 Protein, Human, Recombinant (His)

General Information

Synonyms:	Sodium-dependent vitamin C transporter 2 (hSVCT2);SLC23A1;Nucleobase transporter-like 1 protein;KIAA0238;Yolk sac permease-like molecule 2;SVCT2;Solute carrier family 23 member 2;Na(+)/L-ascorbic acid transporter 2;NBTL1;SLC23A2;YSPL2
Protein Construction:	1-650 aa
Species:	Human
Expression Host:	E. coli
Accession:	Q9UGH3
Molecular Weight:	71.3 kDa (predicted)
AA Sequence:	MMGIGKNTTSKSM EAGSSTEGKYEDEAKHPAFFTLPPVINGGATSSGEQDNEDTELMAIYTTENGIAEKSSLA ETLDSTGSLDPQRSDMIYTIEDVPPWYLCIFLGLQHYLTCSFGTIAVPFLADAMCVGYDQWATSQ LIGTIFFCV GITLLQTTFGCRLPLFQASAFAPARAILS LDKWKCNTT DVS VANGTAELLHTEHIWYPRIREIQGAIIMSSLI EVLIGLLGLPGALLKYIGPLTITPTVALIGLSGFQAAGERAGKHWGIAMLTIFLVLLFSQYARNVKFPLPIYKSKKG WTAYKLQLFKMFPIILAILVSWLLCFIFTVTDVFPDSTKYGFYARTDARQGVLLVAPWFKVPYPFQWGLPTVS AAGVIGMLS AVVASIIESIGDYYACARLSCAPPPPIHAINRGIFVEGLSCVLDGIFGTGNGSTSSSPNIGVLGITKV GSRRVIQCGAALMLALGMIGKFSALFASLPDPVLGALFCTLFGMITAVGLSNLQFIDLN SSRNLFVLGFSIFFGL VLP SYLRQNPLVTGITGIDQVLNVLLTTAMFVGGCVAFILDNTIPGTPEERGIRKWKKG VGGKGNKSLDGMESYN LPFGMNIKKYRCFSYLPISPTFVG YTWKGLRKSDNSRSDSDS QATG

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 90% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in sterile deionized water. The product concentration should not be less than 100 μg/mL. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Sodium/ascorbate cotransporter. Mediates electrogenic uptake of vitamin C, with a stoichiometry of 2 Na(+) for each ascorbate.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481