

## Tetraspanin-7/TSPAN7 Protein, Human, Recombinant (E. coli, His)

### General Information

Synonyms:	Transmembrane 4 superfamily member 2;A15;CD231;Tspan-7;TSPAN7;DXS1692E;Cell surface glycoprotein A15;T-cell acute lymphoblastic leukemia-associated antigen 1 (TALLA-1);Membrane component chromosome X surface marker 1;MXS1;Tetraspanin-7;TM4SF2
Protein Construction:	113-213 aa
Species:	Human
Expression Host:	E. coli
Accession:	P41732
Molecular Weight:	15.6 kDa (predicted)
AA Sequence:	RHEIKDTFLRTYTDAMQTYNGNDERSRAVDHVVQRSLSCCGVQNYTNWSTSPYFLEHGIPPSCCMNETDCNP QDLHNLTVAAATKVNQKGCYDLVTSFMETNM

### QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 90% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	Tris-based buffer, 50% glycerol

### Preparation and Storage

#### Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

#### Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

### Protein Background

May be involved in cell proliferation and cell motility.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481