

Lck Protein, Human, Recombinant (His)

General Information

Synonyms:	Leukocyte C-terminal Src kinase (LSK);LCK;Protein YT16;Tyrosine-protein kinase Lck;T cell-specific protein-tyrosine kinase;Lymphocyte cell-specific protein-tyrosine kinase;Proto-oncogene Lck;p56-LCK
Protein Construction:	1-539 aa
Species:	Human
Expression Host:	E. coli
Accession:	P06239
Molecular Weight:	65.2 kDa (predicted)
AA Sequence:	MGCGCSSHPEDDWMENIDVCENCHYPIVPLDGGKGTLLIRNGSEVRDPLVYEGSNPPASPLQDNLVIALHSYE PSHDGDLGFEKGEQLRILEQSGEWWKAQSLTTGQEGFIPFNFAKANSLEPEPWFFKNLSRKDAERQLLAG NTHGSFLIRESESTAGSFSLSVRDFDQNGQEVVKHYKIRNLDNNGGFYISPRITFPGLHELVRHYTNASDGLCTR LSRPCQTQKPQKPWWEDWEVPRETLKLVRLGAGQFGEVWVGYYNGHTKVAVKSLKQGSMPDAFLAEA NLMKQLQHQLVRLYAVVTQEPIYIITEYMENDTLDSQLEEKGLGASPWGNLQQLLLLPTGSLVDFLKTPS GIKLTINKLLDMAAQIAEGMAFIEERNYIHRDLRAANILVSDTLCKIADFLARLIEDNEYTAREGAKFPIKWTA PEAINYGTFTIKSDVWSFGILLTEIVTHGRIPYPGMTNPEVIQNLERGYRMVRPDNCPEELYQLMRLCWKERPE DRPTFDYLRVLEDFFTATEGQYQPQP

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 90% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	Tris/PBS-based buffer, 6% Trehalose

Preparation and Storage

Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Non-receptor tyrosine-protein kinase that plays an essential role in the selection and maturation of developing T-cells in the thymus and in the function of mature T-cells. Plays a key role in T-cell antigen receptor (TCR)-linked signal transduction pathways. Constitutively associated with the cytoplasmic portions of the CD4 and CD8 surface receptors. Association of the TCR with a peptide antigen-bound MHC complex facilitates the interaction of CD4 and CD8 with MHC class II and class I molecules, respectively, thereby recruiting the associated LCK protein to the vicinity of the TCR/CD3 complex. LCK then phosphorylates tyrosine residues within the immunoreceptor tyrosine-based activation motifs (ITAM) of the cytoplasmic tails of the TCR-gamma chains and CD3 subunits, initiating the TCR/CD3 signaling pathway. Once stimulated, the TCR recruits the tyrosine kinase ZAP70, that becomes phosphorylated and activated by LCK. Following this, a large number of signaling molecules are recruited, ultimately leading to lymphokine production. LCK also contributes to signaling by other receptor molecules. Associates directly with the cytoplasmic tail of CD2, which leads to hyperphosphorylation and activation of LCK. Also plays a role in the IL2 receptor-linked signaling pathway that controls the T-cell proliferative response. Binding of IL2 to its receptor results in increased activity of LCK. Is expressed at all stages of thymocyte development and is required for the regulation of maturation events that are governed by both pre-TCR and mature alpha beta TCR. Phosphorylates other substrates including RUNX3, PTK2B/PYK2, the microtubule-associated protein MAPT, RHOH or TYROBP. Interacts with FYB2.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481