

Measles virus (MeV) (strain Edmonston-Schwarz vaccine) Nucleoprotein/NP Protein (E. coli, His)

General Information

Synonyms: Nucleoprotein;NP;Nucleocapsid protein (NP;Protein N)

Protein Construction: 1-400 aa

Species: MeV

Expression Host: E. coli

Accession: B8PZP3

Molecular Weight: 49.9 kDa (predicted)

AA Sequence:

```
MATLLRSLALFKRNKDKPPITSGSGGAIRGIKHIIIVPIPGDSSITRSRLDRLVRLIGNPDVSGPKLTGALIGILSL
FVESPGQLIQRITDDPDVSIIRLLEVVSQSDQSQSGLTFASRGTNMEDEADQYFSHDDPISSDQSRFGWFGNKEI
SDIEVQDPEGFNMILGTILAQIWWLLAKAVTAPDTAADSELRRWIKYTQRRVVGFEFLERKWLDVVRNRIAE
DLSLRRFMVALILDIKRTPGNKPRIAEMICDIDTYIVEAGLASFILTIKFGIETMYPALGLHEFAGELSTLESMLNLY
QQMGETAPYMVILENSIQNKFSAGSYPLLWSYAMGVGVELENSMGGLNFGRSYFDPAYFRLGQEMVRRSAG
KVSSTLASELGITAEDARLVSEIAMH
```

QC Testing

Biological Activity: Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.

Purity: > 85% as determined by SDS-PAGE.

Endotoxin: < 1.0 EU/μg of the protein as determined by the LAL method.

Formulation: Tris-based buffer, 50% glycerol

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 μg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

A DRUG SCREENING EXPERT

Encapsidates the genome in a ratio of 1 N per 6 ribonucleotides, protecting it from nucleases. The nucleocapsid (NC) has a helical structure with either 12.35 or 11.64 N per turn, approximately 20 nm in diameter, with a hollow central cavity approximately 5 nm in diameter. The encapsidated genomic RNA is termed the NC and serves as template for transcription and replication. During replication, encapsidation by N is coupled to RNA synthesis and all replicative products are resistant to nucleases. N is released in the blood following lysis of measles infected cells, it interacts then with human FCGR2B on immune cells, inducing apoptosis and blocking inflammatory immune response. Ntail binds to a protein on human thymic epithelial cells, termed Nucleoprotein Receptor (NR), inducing growth arrest.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481