

## ADAMTS13 Protein, Mouse, Recombinant (His & Myc)

### General Information

Synonyms:	A disintegrin and metalloproteinase with thrombospondin motifs 13; von Willebrand factor-cleaving protease (vWF-CP; vWF-cleaving protease); ADAMTS-13; ADAM-TS13; Adamts13; Gm710; ADAM-TS 13
Protein Construction:	904-1137 aa
Species:	Mouse
Expression Host:	E. coli
Accession:	Q769J6
Molecular Weight:	32.7 kDa (predicted)
AA Sequence:	WTPLVGLCSISCGRGLKELYFLCMDSVLKMPVQEELCGLASKPPSRWEVCRARPCPARWETQVLAPCPVTCG GGRVPLSVRCVQLDRGHPISVPHSKCSPVPKPGSFEDCSPEPCPARWKVLSLGPSCASCSGLGTATQMVACMQ LDQGHNDNEVNETFCKALVRPQASVPLIADCAFRWHISAWTECSVSCGDGIQRRHDTCLGPQAQVPVPANF CQHLPKPMTVRGCWAGPCA

### QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 85% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	Tris-based buffer, 50% glycerol

### Preparation and Storage

#### Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

#### Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

### Protein Background

Cleaves the vWF multimers in plasma into smaller forms thereby controlling vWF-mediated platelet thrombus formation.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481