

Arginase-1/ARG1 Protein, Mouse, Recombinant (His & SUMO)

General Information

Synonyms: Arg1; Liver-type arginase; Type I arginase; Arginase-1

Protein Construction: 1-323 aa

Species: Mouse

Expression Host: E. coli

Accession: Q61176

Molecular Weight: 50.8 kDa (predicted)

AA Sequence: MSSKPKSLEIIGAPFSKQPRGGVEKGPAALRKAGLLEKLKETEDVDRDHGDLAFVDVPNDSSFQIVKNPRSV
GKANEELAGVVAEVQKNGRVSVVLGGDHLAVGSISGHARVHPDLCVIWVDAHTDINTPLTSSGNLHGQP
VSFLLKELKGFDPVPGFSWVTPCISAKDIVYIGLRDVPGEHYIIKTLGIKYFSMTEVDKLGIGKVMEEFVSYLLG
RKKRPIHLSFDVDGLDPAFTPATGTPVLGGLSYREGLYITEEIKTGLLSGLDIMEVNPTLGKTAEEVKSTVNTAV
ALTACFGTQREGNHKPGTDYLKPPK

QC Testing

Biological Activity: Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.

Purity: > 90% as determined by SDS-PAGE.

Endotoxin: < 1.0 EU/μg of the protein as determined by the LAL method.

Formulation: Tris-based buffer, 50% glycerol

Preparation and Storage

Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Key element of the urea cycle converting L-arginine to urea and L-ornithine, which is further metabolized into metabolites proline and polyamides that drive collagen synthesis and bioenergetic pathways critical for cell proliferation, respectively; the urea cycle takes place primarily in the liver and, to a lesser extent, in the kidneys.;

Functions in L-arginine homeostasis in nonhepatic tissues characterized by the competition between nitric oxide synthase (NOS) and arginase for the available intracellular substrate arginine. Arginine metabolism is a critical regulator of innate and adaptive immune responses. Involved in an antimicrobial effector pathway in polymorphonuclear granulocytes (PMN). Upon PMN cell death is liberated from the phagolysosome and depletes arginine in the microenvironment leading to suppressed T cell and natural killer (NK) cell proliferation and cytokine secretion. In group 2 innate lymphoid cells (ILC2s) promotes acute type 2 inflammation in the lung and is involved in optimal ILC2 proliferation but not survival. Plays a role in the immune response of alternatively activated or M2 macrophages in processes such as wound healing and tissue regeneration, immune defense against multicellular pathogens and parasites, and immune suppression and allergic inflammation; the regulatory outcome seems to be organ specific. In tumor-infiltrating dendritic cells (DCs) and myeloid-derived suppressor cells (MDSCs) plays a role in suppression of T cell-mediated antitumor immunity.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481