

CES1C Protein, Mouse, Recombinant

General Information

Synonyms:	Liver carboxylesterase N;Ces1c;Lung surfactant convertase;Carboxylesterase 1C;PES-N;Es1
Protein Construction:	19-550 aa
Species:	Mouse
Expression Host:	E. coli
Accession:	P23953
Molecular Weight:	58.6 kDa (predicted)
AA Sequence:	<pre>HSELLPPVDDTTQGKVLGKYISLEGFEQPVAVFLGVFPFAKPPLGSLRFAPPQPAEPWSFVKNATSYPPMCSQDA GWAKILSDMFSTEKEILPLKISEDCLYLNIIYSPADLTKSSQLPVMVWIHGGGLVIGGASPYNGLALSAHENVV VTIQYRLGIWGLFSTGDEHSPGNWAHL DQLAALRWVQDNIANFGGNPDSVTIFGESSGGISVSVLVLSPGKD LFHRAISESGVVINTNVGKKNIQAVNEIATLSQCNDTSSAAMVQCLRQKTESELLEISGKLVQYNISLSTMIDG VVLPKAPEEILAEKSFNTVPYIVGFNKQEFGWII PMMLQNLPEGKMNEETASLLRRFHSELNISESMIPAVIE QYLRGVDDPAKKSELILDMFGDIFFGIPAVLLSRSLRDAGVSTYMYEFYRPSFVSDKRPQTVEGDHGD EIFFV FGAPLLKEGASEEETNLSKMVMKFWANFARNGNPNGEGLPHWPEYDEQEGYLQIGATTQQAQRLKAE EEA FWTELLAKNPPETDPTEH</pre>

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 90% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	Tris-based buffer, 50% glycerol

Preparation and Storage

Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

A DRUG SCREENING EXPERT

Involved in the detoxification of xenobiotics and in the activation of ester and amide prodrugs. Involved in the extracellular metabolism of lung surfactant.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481