

TRIM21 Protein, Mouse, Recombinant (His)

General Information

Synonyms:	52 kDa ribonucleoprotein autoantigen Ro/SS-A;Ssa1;52 kDa Ro protein;Tripartite motif-containing protein 21;Ro(SS-A);Ro52;Sjogren syndrome type A antigen (SS-A);Trim21;E3 ubiquitin-protein ligase TRIM21
Protein Construction:	1-470 aa
Species:	Mouse
Expression Host:	Baculovirus Insect Cells
Accession:	Q62191
Molecular Weight:	56.7 kDa (predicted)
AA Sequence:	MSPSTTSKMSLEKMWEEVTCISICLDPMEVPMSEICGHCFCKECIFEVGKNGGSSCPECRQQFLLRNLRPNRHI ANMVENLKQIAQNTKKSTQETHCMKHGKELHLFCEEDGQALCWCAQSGKHRDHTRVPIEEAAKVYQEKIH VVLEKLRKGGKELAEKMEMDLTMQRDQWKRNIQTQKSRIHAEFALQNSLLAQEEQRQLQRLEKQREYLRLLG KKEAELAENKQALQELISELERRIRGSELELLQEVRIILERSGSWNLDLTDIDAPDLTSTCPVPRKKMLRWCW HITLDRNTANSWLIISKDRRQVRMGDTHQNVSDNKERFSNYPMVLGAQRFSSGKMYWEVDVTQKEAWDLG VCRDSVQRKGQFSLSPENGFWTIWLWQDSYEAGTSPQTTLHIQVPPCQIGIFVDYEAGVVSFYNITDHGSLIYT FSECVFAGPLRPFNVGFNYSGGNAAPLKLCLPKM

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 85% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in sterile deionized water. The product concentration should not be less than 100 μg/mL. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

E3 ubiquitin-protein ligase whose activity is dependent on E2 enzymes, UBE2D1, UBE2D2, UBE2E1 and UBE2E2. Forms a ubiquitin ligase complex in cooperation with the E2 UBE2D2 that is used not only for the ubiquitination of USP4 and IKBKB but also for its self-ubiquitination. Component of cullin-RING-based SCF (SKP1-CUL1-F-box protein) E3 ubiquitin-protein ligase complexes such as SCF(SKP2)-like complexes. A TRIM21-containing SCF(SKP2)-like complex is shown to mediate ubiquitination of CDKN1B ('Thr-187' phosphorylated-form), thereby promoting its degradation by the proteasome. Monoubiquitinates IKBKB that will negatively regulates Tax-induced NF-kappa-B signaling. Negatively regulates IFN-beta production post-pathogen recognition by polyubiquitin-mediated degradation of IRF3. Mediates the ubiquitin-mediated proteasomal degradation of IgG1 heavy chain, which is linked to the VCP-mediated ER-associated degradation (ERAD) pathway. Promotes IRF8 ubiquitination, which enhanced the ability of IRF8 to stimulate cytokine genes transcription in macrophages. Plays a role in the regulation of the cell cycle progression. Enhances the decapping activity of DCP2. Exists as a ribonucleoprotein particle present in all mammalian cells studied and composed of a single polypeptide and one of four small RNA molecules. At least two isoforms are present in nucleated and red blood cells, and tissue specific differences in RO/SSA proteins have been identified. The common feature of these proteins is their ability to bind HY RNAs.2. Involved in the regulation of innate immunity and the inflammatory response in response to IFNG/IFN-gamma. Organizes autophagic machinery by serving as a platform for the assembly of ULK1, Beclin 1/BECN1 and ATG8 family members and recognizes specific autophagy targets, thus coordinating target recognition with assembly of the autophagic apparatus and initiation of autophagy. Acts as an autophagy receptor for the degradation of IRF3, hence attenuating type I interferon (IFN)-dependent immune responses. Represses the innate antiviral response by facilitating the formation of the NMI-IFI35 complex through 'Lys-63'-linked ubiquitination of NMI.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481