

NEFL Protein, Mouse, Recombinant (His)

General Information

Synonyms:	Nfl;Nefl;Nf68;68 kDa neurofilament protein;Neurofilament light polypeptide;NF-L; Neurofilament triplet L protein
Protein Construction:	2-543 aa
Species:	Mouse
Expression Host:	P. pastoris (Yeast)
Accession:	P08551
Molecular Weight:	63.4 kDa (predicted)
AA Sequence:	<pre>SSFGYDPPYFSTSYKRRYVETPRVHISSVRSYGSTARSAYSSYSAPVSSLSVRRSYSSSSGSLMPLENLDLSQV AAISNDLKSIRTQEKAQLQDLNDRFASFIERVHELEQQNKVLEAELLVLRQKHSEPSRFRALYEQEIIDLRLAA EDATNEKQALQGEREGLEETLRNLQARYEEEEVLSREDAEGRLMEARKGADEAALARAEELEKRIDSLMDEIAFL KKVHEEEI AELQAQIQYAQISVEMDVSSKPDLSAALKDIRAQYEKLAAKNMQNAEEWFKSRFTVLTESAANKT DAVRAAKDEVSESRLLKAKTLEIEACRGMNEALEKQLQELEDKQNADISAMQDTINKLENELRSTKSEMARY LKEYQDLLNVKMALDIEIAAYRKLLERGEETRLSFTSVGSITSGYSQSSQVFGRSAYSLQSSSYLMSARSFPAYY TSHVQEEQTEVEETIEATKAEAAKDEPPSEGEAEKEEKEEKEEKEEKEEKEEKEEKEEKEEKEEKEEKEEKEE DTKESEEEKEEKEESAGEEQVAKKKD</pre>

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 90% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	Tris-based buffer, 50% glycerol

Preparation and Storage

Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Neurofilaments usually contain three intermediate filament proteins: NEFL, NEFM, and NEFH which are involved in the maintenance of neuronal caliber. May additionally cooperate with the neuronal intermediate filament proteins PRPH and INA to form neuronal filamentous networks.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481