

Nidogen-1 Protein, Mouse, Recombinant (E. coli, His)

General Information

Synonyms: Ent;Nid1;NID-1;Nidogen-1;Entactin

Protein Construction: 428-665 aa

Species: Mouse

Expression Host: E. coli

Accession: P10493

Molecular Weight: 30.3 kDa (predicted)

AA Sequence: SPQRVNGKVKGRIFVGSSQVPVVFENTDLHSYVVMNHGRSYTAISTIPETVGYSLPLAPIGGIIGWMFAVEQD
GFKNGFSITGGEFTRQAEVTFLGHPGKLVKQQFSGIDEHGHLTISTELEGRVPQIPYGASVHIEPYTELYHYSS
SVITSSSTREYTMPEPDQDGAAPSHTHIYQWRQTITFQECAHDDARPALPSTQQLSVDSVFLYNKEERILRY
ALSNSIGPVRDGPDA

QC Testing

Biological Activity: Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.

Purity: > 90% as determined by SDS-PAGE.

Endotoxin: < 1.0 EU/μg of the protein as determined by the LAL method.

Formulation: Tris-based buffer, 50% glycerol

Preparation and Storage

Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Sulfated glycoprotein widely distributed in basement membranes and tightly associated with laminin. Also binds to collagen IV and perlecan. It probably has a role in cell-extracellular matrix interactions.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481