

Jagged 1 Protein, Mouse, Recombinant (His & Myc & SUMO)

General Information

| | |
|-----------------------|--|
| Synonyms: | Protein jagged-1;CD339;Jagged1;Jag1 |
| Protein Construction: | 33-334 aa |
| Species: | Mouse |
| Expression Host: | E. coli |
| Accession: | Q9QXX0 |
| Molecular Weight: | 53.6 kDa (predicted) |
| AA Sequence: | GQFELEILSMQNVNDELQNGNCCGGVRNPGDRKCTRDECDTYFKVCLKEYQSRVTAGGPCSFSGSGSTPVIGG NTFNLKASRGNDRNRIVLPFSFAWPRSYYLLVEAWDSSNDTIQPDSIIKASHSGMINPSRQWQTLKQNTGIA HFYQIRVTCDDHYYGFGCNKFCRPRDDFFGHYACDQNGNKTCMEGWMGPDCNKAICRQGCSPKHGSKL PGDCRCQYGWQGLYCDKCIHPGCVHGTCEPQWQCLCETNWGGQLCDKDLNYCGTHQCPLNRGTCSNTG PDKYQCSCPEGYS GPNCE |

QC Testing

| | |
|----------------------|---|
| Biological Activity: | Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first. |
| Purity: | > 90% as determined by SDS-PAGE. |
| Endotoxin: | < 1.0 EU/μg of the protein as determined by the LAL method. |
| Formulation: | Tris-based buffer, 50% glycerol |

Preparation and Storage

Reconstitution:
Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 μg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Ligand for multiple Notch receptors and involved in the mediation of Notch signaling. May be involved in cell-fate decisions during hematopoiesis. Seems to be involved in early and late stages of mammalian cardiovascular

development. Inhibits myoblast differentiation. May regulate fibroblast growth factor-induced angiogenesis.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481