

## Thrombospondin-2/THBS2 Protein, Mouse, Recombinant (E. coli, His)

### General Information

Synonyms:	Thrombospondin-2;Tsp2;Thbs2
Protein Construction:	19-232 aa
Species:	Mouse
Expression Host:	E. coli
Accession:	Q03350
Molecular Weight:	28.1 kDa (predicted)
AA Sequence:	GDHVKDTSFDFLSISNINRKTIGAKQFRGPDGPVPAYRFVRFDYIPPVNTDDLNRIVKLARRKEGFFLTAQLKQDRKSRGTLVLEGPQTSQRQFEIVSNGPGDTLDLNYWVEGNQHTNFLEDVGLADSQWKNTVQVASDTYSLYVGCDLIDSVTLEEPFYEQLEVDRSRMYVAKGASRESHFRGLLQNVHLVFADSVEDILSKKGCQHSQGA

### QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 90% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	Tris-based buffer, 50% glycerol

### Preparation and Storage

**Reconstitution:**  
A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

**Stability & Storage:**

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

**Shipping:**

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

### Protein Background

Adhesive glycoprotein that mediates cell-to-cell and cell-to-matrix interactions. Ligand for CD36 mediating antiangiogenic properties.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481