

## TNNI1 Protein, Mouse, Recombinant (Flag & His & Myc)

### General Information

Synonyms:	Troponin I, slow skeletal muscle;Troponin I, slow-twitch isoform;slow skeletal muscle;Tnni1
Protein Construction:	2-187 aa
Species:	Mouse
Expression Host:	E. coli
Accession:	Q9WUZ5
Molecular Weight:	29.6 kDa (predicted)
AA Sequence:	PEVERKSKITASRKLMLKSLMLAKAKECWEQEHEEREAEKVRYLSERIPTLQTRGLSLSALQDLCRELHAKVEV VDEERYDIEAKCLHNTREIKDLKLVLDLRGKFKRPPLRRVRVSADAMLRALLGSKHKVSMDLRANLKSVKKE DTEKERPVEVGDWRKNVEAMSGMEGRKKMFDAAKSPTSQ

### QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 85% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	Tris-based buffer, 50% glycerol

### Preparation and Storage

#### Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

#### Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

### Protein Background

Troponin I is the inhibitory subunit of troponin, the thin filament regulatory complex which confers calcium-sensitivity to striated muscle actomyosin ATPase activity.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286   E\_mail: info@targetmol.com   Address: 34 Washington Street, Wellesley Hills, MA 02481