

Galectin-4/LGALS4 Protein, Rat, Recombinant (His & SUMO)

General Information

Synonyms:	Galectin-4;L-36 lactose-binding protein (L36LBP);Gal-4;Lgals4;Lactose-binding lectin 4
Protein Construction:	1-324 aa
Species:	Rat
Expression Host:	E. coli
Accession:	P38552
Molecular Weight:	52.3 kDa (predicted)
AA Sequence:	MAYVPAPGYQPTYNPTLPYKRPIPGGLSVGMSIYIQGIKDNMRRFHVNFVAVGQDEGADIAFHFNPRFDGW DKVVFNTMQSGQWGKEEKKKSMFPQKGGHFFELVFMVMSEHYKVVVNGTPFYEYGHRLPLQMVTHLQVDG DLELQSINFLGGQPAASQYPGTMTIPAYPSAGYNPPQMNSLPVMAGPPIFNPPVPYVGTLQGGLTARRTIIKKG YVLPTAKNLIINFKVGSTGDIAFHMNPRIGDCVVRNSYMNGSWGSEERKIPYNPFGAGQFFDLSIRCGTDRFK VFANGQHLLDFSHRFQAFQRVDMLEIKGDITLSYVQI

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 90% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	Tris-based buffer, 50% glycerol

Preparation and Storage

Reconstitution:
A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Galectin that binds lactose and a related range of sugars.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481