

## Serpin A3K Protein, Rat, Recombinant (His &amp; Myc)

## General Information

Synonyms:	Spin2b;Contrapsin-like protease inhibitor 1;GHR-P63;CPI-21;Serpina3k;Kallikrein-binding protein (KBP);Serine protease inhibitor 2 (SPI-2);SPI-2.3;Serpin A3K;Thyroid hormone-regulated protein;Growth hormone-regulated proteinase inhibitor;Serine protease inhibitor A3K
Protein Construction:	21-416 aa
Species:	Rat
Expression Host:	E. coli
Accession:	P05545
Molecular Weight:	52.0 kDa (predicted)
AA Sequence:	DGILGRDTLPHEHQGKGRQLHSLTLASINTDFTLSLYKKLALRNPDKNVVFSPLSISAALAILSLGAKDSTMEIILEVKFNLTETEEEEIHQFGHLLQRLSQPEDQAEINTGSALFIDKEQPILSEFQEKTRALYQAEAFVADFKQCNEAKKFINDYVSNQTQGKIAELFSELDERTSMVLVNYLLFKGKWKVPFNPNDTFESEFYLDEKRSVKVPMMKIKDLTPYIRDEELSCSVLELKYTGNASALFILPDQGKMQQVESSLQPETLKKWKDSLPRRIISELRMPKFSISTDYNLEEVLPGLGIRKIFSQQADLSRITGTKNLHVSQVVKAVLDVDETGTGEGAAATAVTAALKSLPQTIPLLNFRPFMLVITDNNQSVFFMGKVTNPM

## QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 85% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/µg of the protein as determined by the LAL method.
Formulation:	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.

## Preparation and Storage

## Reconstitution:

Reconstitute the lyophilized protein in sterile deionized water. The product concentration should not be less than 100 µg/mL. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

## Stability &amp; Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

## Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Binds to and inhibits kallikreins. Inhibits trypsin but not chymotrypsin or elastase.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481