

## Rotavirus A (strain RVA/Cow/France/RF/1975/G6P6[1]) VP6 Protein (His)

### General Information

Synonyms: Intermediate capsid protein VP6

Protein Construction: 1-397 aa

Species: RV-A

Expression Host: E. coli

Accession: P04509

Molecular Weight: 48.8 kDa (predicted)

AA Sequence:

```
MDVLYSLSKTLKDARDKIVEGTLYSNVSDLIQQFNQMIITMNGNEFQTGGIGNLPIRNWNFDFGLLGTTLLNLD
ANYVETARNTIDYFVDFVDNVCMDVMVRESQRNGIAPQSDSLIKLSGIKFKRINFDNSSEYIENWNLQNRQR
TGFTFHKPNIFPYSASFTLNRSQPAHDNLMGTMWLNAGSEIQVAGFDYSCAINAPANTQQFEHIVQLRRVLT
TATITLLPDAERFSFPRVITSADGATTWYFNPVILRPNNVEIEFLLNGQIINTYQARFGTIIARNFDTIRLSFQLMR
PPNMTPAVAALFPNAQPFEHHATVGLTLRIESAVCESVLADASETMLANVTSVRQEYAI PVGPFVPPGMNWT
DLITNYSRENDLQRVFTVASIRSMVLK
```

### QC Testing

Biological Activity: Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.

Purity: > 90% as determined by SDS-PAGE.

Endotoxin: < 1.0 EU/μg of the protein as determined by the LAL method.

Formulation: Tris-based buffer, 50% glycerol

### Preparation and Storage

Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months.

Please avoid multiple freeze-thaw cycles and store products in aliquots.

*Actual storage temperature shall be subject to the COA.*

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

### Protein Background

Intermediate capsid protein that self assembles to form an icosahedral capsid with a T=13 symmetry, which consists of 230 trimers of VP6, with channels at each of its five-fold vertices. This capsid constitutes the middle

concentric layer of the viral mature particle. The innermost VP2 capsid and the intermediate VP6 capsid remain intact following cell entry to protect the dsRNA from degradation and to prevent unfavorable antiviral responses in the host cell during all the replication cycle of the virus. Nascent transcripts are transcribed within the structural confines of this double-layered particle (DLP) and are extruded through the channels at the five-fold axes. VP6 is required for the transcription activity of the DLP.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481