

CAK1 Protein, *S. cerevisiae*, Recombinant (E. coli, His & Myc)

## General Information

Synonyms:	CDK-activating kinase; CIV1; CAK1; Serine/threonine-protein kinase CAK1
Protein Construction:	1-368 aa
Species:	<i>Saccharomyces cerevisiae</i>
Expression Host:	<i>E. coli</i>
Accession:	P43568
Molecular Weight:	49.6 kDa (predicted)
AA Sequence:	MKLDSIDITHCQLVKSTRTARIYRSDTYAIKCLALDFDIPPHNAKFEVSILNKLGNKCKHILPLLESKATDNNDDL LLFPFEEMNLYEFMQMHYKRDRRKNPYYDLLNPSIPIVADPPVQKYTNQLDVNRYSLSFRRQMVEGIAFLHE NKIIHRDIKPQNIMLTNNTSTVSPKLYIIDFGISYDMANNSQTSAPMDSKVTDISTGIYKAPEVLFGVKCYDGG VDVWSLLIISQWFQRETSRMGHVPAMIDDGSDDMNSDGSDFRLICSI FEKLGIPSIQKWEEVAQHGSVDAFV GMFGADGDGKYVLDQEKDQVQISIVERNMPRLDEIADV KVKQKFINCILGMVSFSPNERWSCQRILQELEKP

## QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 85% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	Tris-based buffer, 50% glycerol

## Preparation and Storage

## Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

## Stability &amp; Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

## Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

**This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use**

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481