

## Varicella-zoster virus (strain Dumas) glycoprotein E/gE Protein (His &amp; SUMO)

## General Information

Synonyms: gE;Envelope glycoprotein E

Protein Construction: 31-623 aa

Species: HHV-3

Expression Host: E. coli

Accession: P09259

Molecular Weight: 85.4 kDa (predicted)

AA Sequence: SVLRYDDFHTDEDKLDTNSVYEPYYHSDHAESSWVNRGESSRKAYDHNSPYIWPRNDYDGFLENAHEHHG  
VYNQGRGIDSGERLMQPTQMSAQEDLGDGTGIHVPTLNGDDRHKIVNVDQRQYGDVFKGDLNPKPQGQR  
LIEVSVEENHPFTLRAPIQRIYGVRYTETWSFLPSLTCTGDAAPAIQHICLKHTTCFQDVVVDVCAENTKEDQL  
AEISYRFQGGKEADQPWIVVNTSTLFDELELDPPEIEPGVLKVLRTKQYLGVIWVNMVNRGSDGTSTYATFLVTW  
KGDEKTRNPTPAVTPQPRGAEFHMWNYHSHVFSVGDTFSLAMHLQYKIHEAPFDLLEWLYVPIDPTCQPM  
RLYSTCLYHPNAPQCLSHMNSGCTFTSPHLAQRVASTVYQNCHEADNYTAYCLGISHMEPSFGLLHDGGTTL  
KFVDTPESLSGLYVFVVFNGHVEAVYTVVSTVDHFVNAIEERGFPTAGQPPATTKPKEITPVNPGTSPLLR  
YAAWTGGLAAVVLLCLVIFLICTAKRMRVKAYRVDKSPYNQSMYYAGLPVDDFEDSESTDEEEFGNAIGGSH  
GGSSYTVYIDKTR

## QC Testing

Biological Activity: Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.

Purity: > 85% as determined by SDS-PAGE.

Endotoxin: < 1.0 EU/μg of the protein as determined by the LAL method.

Formulation: Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0

## Preparation and Storage

Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months.

Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

### Protein Background

Envelope glycoprotein that binds to the potential host cell entry receptor IDE.; In epithelial cells, the heterodimer gE/gI is required for the cell-to-cell spread of the virus, by sorting nascent virions to cell junctions. Once the virus reaches the cell junctions, virus particles can spread to adjacent cells extremely rapidly through interactions with cellular receptors that accumulate at these junctions. Implicated in basolateral spread in polarized cells. In neuronal cells, gE/gI is essential for the anterograde spread of the infection throughout the host nervous system. Together with US9, the heterodimer gE/gI is involved in the sorting and transport of viral structural components toward axon tips.; The heterodimer gE/gI serves as a receptor for the Fc part of host IgG. Dissociation of gE/gI from IgG occurs at acidic pH. May thus be involved in anti-VZV antibodies bipolar bridging, followed by intracellular endocytosis and degradation, thereby interfering with host IgG-mediated immune responses.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481